HemoCue® Hb 301

Operating manual

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Manual de funcionamiento

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GB/US

Table of Contents

HemoCue Hb 301 system.....	5
Components.....	6
Start up.....	8
Set-up.....	10
Audio signal.....	10
Measuring.....	14
Capillary blood.....	14
Measuring.....	22
Control material, venous or arterial blood.....	22
Set-up.....	24
PC/Printer function.....	24
Maintenance.....	26
Troubleshooting Guide.....	28
Specifications.....	44

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Sistema HemoCue® Hb 301

GB/US

Thank you for choosing the HemoCue® Hb 301 system. The system is designed for quantitative point-of-care whole blood hemoglobin determination in primary care or blood donation settings using a specially designed analyzer, the HemoCue Hb 301 Analyzer, and specially designed microcuvettes, the HemoCue Hb 301 Microcuvettes. The HemoCue Hb 301 system is for *In Vitro* Diagnostic use only. Capillary, venous or arterial whole blood may be used. The test results can be printed or down-loaded to a PC.

All system components are designed and manufactured to provide maximum safety. Any other use of the system may impair the safety.

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Components
Composants
Componentes
Componentes

GB/US

1. HemoCue Hb 301 Analyzer*
2. AC adapter**
3. 4 type AA batteries***
4. A vial of HemoCue Hb 301 Microcuvettes
5. HemoCue Hb 301 Operating Manual, HemoCue Hb 301 Quick Reference Guide and Instruction CD
6. HemoCue Cleaner

Open the carton and lift out the instrument and accessories.

If no AC power is available, use 4 type AA batteries. On the bottom of the analyzer there is a lid covering the battery compartment. Press the flap to remove the lid. Place the batteries in the battery compartment and replace the lid.

Consult local environmental authorities for proper disposal of batteries.

* Do not open the cover of the analyser.

**  Only use adapters listed under specifications.

*** Not included.

For information about HemoCue Hb 301 Microcuvettes please contact your HemoCue distributor.

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
Funcionamento

GB/US

Only use adapters approved by HemoCue, as listed under "Technical Specifications".

1. If AC power is available, use the power adapter.
2. If no power is available, insert four type AA batteries, 1.5 V.
3. Pull the cuvette holder out to the loading position.
4. Press and hold the button until the display is activated (all symbols appear on the display). The analyzer performs a self-test, and after approximately 10 seconds the display will show three flashing dashes and the HemoCue symbol. This indicates that the analyzer is ready for use.
5. If the cuvette holder is in the loading position, the most recent result will be displayed.
If the power adapter had previously been disconnected for more than 5 minutes, and there are no batteries in the analyzer, "Hb" will be displayed instead of the most recent result.

To turn the analyzer off, press and hold the button until the display reads OFF and then goes blank. When operating on battery power, the analyzer will automatically turn off after approximately 5 minutes.

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
Sinal sonoro

GB/US

If the audible signal is activated, a signal will be heard when the measurement is completed, or if an error code is displayed.

1. Make sure that the analyzer is turned off. Press the On/Off button for approximately 10 seconds.
2. The display now shows a flashing bell-symbol.
3. Press the button rapidly to change between signal on and off.

When the setting is completed, press the button for approximately 3 seconds until the bell stops flashing. The analyzer will return to the ready mode.

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GB/US

In case of peripheral circulatory failure, hemoglobin measurements from capillary samples may be misleading. Confirmation of an unacceptable or unexpected laboratory result is considered as a common, well established routine in the laboratory. Extra cellular fluid present in capillary samples due to individual physiology or the affects of unacceptable capillary sampling technique can cause false hemoglobin results.

One of several pre analytical factors associated with technique revolves around the size and action mode of the lancet and how it is used. This and other factors associated with capillary sampling technique can have an impact on the puncture, blood flow and ultimately on the results obtained. Additional factors unrelated to sampling technique, such as an individual's capillary circulation status, must also be recognized as potential causes of false results. Therefore when testing for hemoglobin using capillary blood as the sample, confirmation of an unacceptable or unexpected result is needed to rule out potential pre analytical factors as the cause. This is particularly helpful when a cut off value is used, i.e. when screening for nutritional status based on hemoglobin. Since the individual performing the test may not always be aware of all the potential pre analytical factors, or even those that are directly or indirectly associated with sampling techniques, a second finger stick should be performed when wanting to rule out pre analytical factors as a cause of an unacceptable or unexpected result. Performance of the second finger stick by a different individual may aid in lessening the affect of the factors that could have been technique related.

Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.

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
GB/US

Note. Always handle blood specimens with care, as they might be infectious. Always wear protective gloves when handling blood specimens.

1. To perform a test, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol. Remove a HemoCue Hb 301 Microcuvette from the vial.
2. Make sure the patient's/blood donor's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean with disinfectant and allow to dry completely or wipe off with a dry, lint-free wipe.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip.

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


GB/US

5. Sample at the side of the fingertip.
6. While applying light pressure towards the fingertip, puncture the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears.

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


11

GB/US

9. When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill! If a second sample is to be taken, it is important that this be done after the measurement of the first sample is complete.
10. Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe. Do not touch the open end of the microcuvette.
11. Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a new drop of blood. Small bubbles around the edge can be ignored.

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



16

GB/US

12. Place the filled microcuvette in the cuvette holder within 40 seconds after filling.
13. Gently touch the cuvette holder. It will automatically slide to the measuring position.
14. During the measurement “” will be shown.
15. After approximately 10 seconds, the hemoglobin value is displayed. The result will remain on the display as long as the cuvette holder is in the measuring position. Do not remeasure the cuvette. When operating on battery power, the analyzer will automatically turn off after approximately 5 minutes.
16. Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal.

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
Measuring Control material, venous or arterial blood

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GB/US

Note. Always handle blood specimens with care, as they might be infectious. Always wear protective gloves when handling blood specimens.

1. If the specimen has been stored in a refrigerator, allow it to warm up to room temperature before mixing. The blood or control material should be well mixed prior to measurement.
2. Place a drop of blood or control material onto a hydrophobic plastic or glass slide using a pipette or other suitable transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill. Wipe off excess sample from the outside of the microcuvette with a clean, lint-free wipe. Do not touch the open end of the microcuvette. If a second sample is to be taken, it is important that this be done after the measurement of the first sample is complete.
4. Look for air bubbles in the filled microcuvette. If present, discard it and fill a new microcuvette from a new drop of blood or control sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–16 on page 21.

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
Maintenance
Maintenance
Mantenimiento
Manutenção

GB/US

The cuvette holder should be cleaned after each day of use. Make sure that the analyzer is turned off and that the display is blank.

1. Pull the cuvette holder out to its loading position. Carefully press the small catch on the cuvette holder using a pointed object.
2. While pressing the catch, carefully pull out the cuvette holder.
3. Clean the cuvette holder with alcohol (70 %) or mild detergent.
4. If the optronic unit becomes dirty, an error code will be displayed. To clean the optronic unit, push a HemoCue Cleaner* swab into the opening of the cuvette holder. Move from side to side 5–10 times. If the swab is stained, repeat with a new swab. Clean the black details of the cuvette holder housing (4a). Wait 15 minutes before replacing the cuvette holder. Make sure that the cuvette holder is locked into place by the small catch. The cover may be cleaned with alcohol (70 %) or a mild detergent.

*One HemoCue Cleaner is included with the analyzer. Contact your local distributor to order more.

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Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact your local HemoCue distributor or HemoCue AB. The analyzer should be cleaned as recommended under "Maintenance" prior to service or disposal. Consult local environmental authorities for proper disposal. The analyzer has no serviceable parts.

Symptom	Explanation	Action
The analyzer shows an error code.	May be a temporary fault.	Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
E00	No stable endpoint of the measurement is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact your distributor.
E01–E05	1. Dirty optronic unit or faulty electronics or optronic unit.	1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact your distributor.
E06	1. Unstable blank value. The analyzer might be cold.	1. Turn off the analyzer and allow it to reach room temperature. If the problem continues, the analyzer needs service. Contact your distributor.
E07	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries, 4 type AA. 1b. Use the power adapter.
E08	The absorbance is too high. 1. Light blocking item in the cuvette holder	1a. Check that the analyzer and microcuvettes are used according to the HemoCue Hb 301 operating manual and instructions for use. 1b. The analyzer needs service. Contact your distributor.

Symptom	Explanation	Action
E10–E30	<ol style="list-style-type: none"> 1. Dirty optronic unit or faulty electronics or optronic unit. 	<ol style="list-style-type: none"> 1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact your distributor.
E40	<ol style="list-style-type: none"> 1. The cuvette holder is not replaced properly after cleaning. 2. Dirty optronic unit. 3. The microcuvette is not a HemoCue Hb 301 Microcuvette. 4. The microcuvette is damaged. 	<ol style="list-style-type: none"> 1. Make sure that the cuvette holder is replaced properly as described in the maintenance section. 2. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 3. Only use HemoCue Hb 301 Microcuvettes in the HemoCue Hb 301 Analyzer. 4. Take a new microcuvette and repeat the measurement.
E41	<ol style="list-style-type: none"> 1. Incorrect adapter. 2. Bad batteries. 3. Hardware error. 	<ol style="list-style-type: none"> 1. Use an adapter listed under Adapters in Specifications. 2. Replace batteries. 3. The analyzer needs service. Contact your distributor.
E42–49	<ol style="list-style-type: none"> 1. The optronic unit has been scratched due to incorrect maintenance. 2. Hardware error. 	<ol style="list-style-type: none"> 1. Clean the optronic unit as described in the maintenance section, using the HemoCue Cleaner. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor.
HHH	<ol style="list-style-type: none"> 1. Measured value exceeds 25.6 g/dL (256 g/L, 15.9 mmol/L). 	
No characters on the display.	<ol style="list-style-type: none"> 1. The analyzer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. 	<ol style="list-style-type: none"> 1a. Check that the power adapter is connected to the analyzer and the AC power supply. 1b. Check that the cable is not damaged. 2. Turn off the analyzer and replace the batteries, 4 type AA. 3. The analyzer needs service. Contact your distributor.
The display contains erroneous characters.	<ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor.

Symptom	Explanation	Action
The display shows "FIR".	This function is for manufacturing use only.	<ol style="list-style-type: none"> 1a. Remove and replace all cables and/or batteries, and restart. 1b. The analyzer needs service. Contact your distributor.
The display shows "▣▣▣".	<ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on AC power, the power adapter or the circuit board is out of order. 	<ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, 4 type AA. 2a. Check that the power adapter is properly connected and working. 2b. The analyzer needs service. Contact your distributor.
The display does not switch from "Σ" and "Hb" to three flashing dashes and "●" (ready for measuring).	<ol style="list-style-type: none"> 1. The cuvette holder sensor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor.
Measurement on control materials are out of range- either too high or too low.	<ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The controls are beyond their expiration dates or have been improperly stored. 4. The control has not been mixed properly and/or is not at room temperature. 5. The microcuvette has not been placed in the analyzer within 40 seconds of filling. 6. Air bubbles in the microcuvette. 7. The optronic unit is dirty. 8. The control is not suitable for use with the HemoCue Hb 301 system. 9. The calibration of the analyzer has been changed. 	<ol style="list-style-type: none"> 1. Check the expiration date and the storage conditions of the microcuvettes. 2. Remeasure the control with a new microcuvette. 3. Check the expiration date and the storage conditions of the control. Remeasure the control with a new microcuvette. If the problem continues , contact the manufacturer of the control. 4. Make sure that the control is mixed properly and at room temperature. If the problem continues , contact the manufacturer of the control. 5. Remeasure the control with a new microcuvette. 6. Check the microcuvette for air bubbles. Remeasure the control with a new microcuvette. 7. Clean the optronic unit as described in the maintance section. 8. Contact your distributor for control information. 9. The analyzer needs service. Contact your distributor.

Symptom	Explanation	Action
<p>Measurement on patient samples are higher or lower than anticipated.</p>	<ol style="list-style-type: none"> 1. Improper sampling technique. 2. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 3. The optical eye of the microcuvette is contaminated. 4. Air bubbles in the microcuvette. 5. The optronic unit is dirty. 6. The calibration of the analyzer has been changed. 	<ol style="list-style-type: none"> 1. See pages 10–19 in this manual. 2. Check the expiration date and the storage conditions of the microcuvettes. 3. Remeasure the sample with a new microcuvette. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit as described in the maintenance section. 6. The analyzer needs service. Contact your distributor.

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Specifications

Intended Purpose/Intended Use

The HemoCue® Hb 301 system is designed for quantitative point-of-care hemoglobin determination in primary care or blood donation settings using a specially designed analyzer, the HemoCue Hb 301 Analyzer, and specially designed microcuvettes, the HemoCue Hb 301 Microcuvettes. The HemoCue Hb 301 system is for *In Vitro* Diagnostic use only. The HemoCue Hb 301 Analyzer is only to be used with HemoCue Hb 301 Microcuvettes.

Principle of the method/procedure

The system consists of an analyzer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette. A blood sample of approximately 10 µL is drawn into the cavity by capillary action. The measurement takes place in the analyzer, which measures the absorbance of whole blood at a Hb/HbO₂ isobestic point. The analyzer measures at two wavelengths (506 and 880 nm) in order to compensate for turbidity. The HemoCue Hb 301 System is calibrated against the hemoglobincyanide (HiCN) method, the international reference method for the determination of the hemoglobin concentration in blood³. The system is factory calibrated and needs no further calibration.

Composition

The microcuvette is made of polystyrene plastic and contains no active ingredients.

Storage and handling

HemoCue® Hb 301 Microcuvettes

The microcuvettes are to be stored at 10–40 °C (50–104 °F). Once the seal of the vial is broken, the microcuvettes are stable for 3 months. An unopened vial of microcuvettes can be stored for a shorter period of time (6 weeks) between -18–50 °C (0.4–122 °F). Use the microcuvettes prior to the expiration date that is printed on the package. Keep the vial properly closed. All unused microcuvettes should remain in the original package.

HemoCue® Hb 301 Analyzer

The analyzer can be stored and transported in temperatures between 0–50 °C (32–122 °F). Operating temperature is 10–40 °C (40–104 °F). Allow the analyzer to reach ambient temperature before use. The analyzer should not be stored or operated at < 5 % or > 90 % non-condensing humidity.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. Appropriate anticoagulant (eg EDTA or heparin) may be used, preferably in solid form to avoid dilutional effects. Tubes containing fluoride should not be used. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing².

Materials required

- HemoCue Hb 301 Analyzer
- HemoCue Hb 301 Microcuvettes
- Lancet (capillary samples)
- Pipette or other transfer device (venous or arterial samples)
- Lint-free tissue (non-fraying)

Quality Control

The HemoCue Hb 301 Analyzer has an internal quality control the "selftest". Every time the analyzer is turned on, it will automatically verify the measurement performance. This test is performed at regular intervals if the analyzer remains switched on. Upon passing the selftest, the display will show the HemoCue symbol and three flashing dashes, indicating that the analyzer is ready to perform a measurement. An error code will be displayed if the selftest fails. If quality control checks are required for regulatory reasons, contact HemoCue or your distributor for current recommendations for controls. Please refer to local guidelines for recommended frequency of use.

Expected values²

Children 11.0–14.0 g/dL (110–140 g/L, 6.8–8.7 mmol/L)
Woman 12.0–15.0 g/dL (120–150 g/L, 7.4–9.3 mmol/L)
Men 13.0–17.0 g/dL (130–170 g/L, 8.1–10.5 mmol/L)
Children, 2 years to teenage, gradually increase to adult values.
Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Measuring range

Displayed range: 0–25.6 g/dL (0–256 g/L, 0–15.9 mmol/L). The method used in the HemoCue Hb 301 system has according to the NCCLS document EP6 been demonstrated to be linear between 2–25 g/dL, (within 4 % difference in 2–5 g/dL and within 3 % difference in 5–25 g/dL)⁵.

Limitations of the procedure

Refer to the HemoCue Hb 301 Microcuvettes package insert¹.

Specific performance characteristics

Within-run and Total precision

Within-run and total precision was determined according to the NCCLS Document EP5-A⁴. The results given below in “Within-run precision” and “Total precision” come from 3 batches of HemoCue Hb 301 Microcuvettes, 5 HemoCue Hb 301 Analyzers and 5 operators. No recalibration was performed during the analyzing period. Commercially available controls at 3 different levels were used. The hemoglobin concentration was measured in duplicate twice a day, morning and afternoon, during 20 consecutive days.

Control level	N	x ^x g/dL	Within-run Precision		Total Precision	
			SD g/dL	CV %	SD g/dL	CV %
1	400	7.3	0.059	0.82	0.066	0.91
2	400	13.2	0.106	0.80	0.122	0.92
3	400	17.2	0.135	0.78	1.152	0.88

Accuracy

The results of the comparison studies between the HemoCue Hb 301 system and the International Council for Standardization in Haematology method (ICSH)³ are summarized in the table below. The study was performed on 3 batches of HemoCue Hb 301 Microcuvettes, 4 HemoCue Hb 301 Analyzers and by 8 operators. No recalibration was performed during the analyzing period.

Study	N	Min g/dL	Max g/dL	Regression line	Correlation coefficient (r)
1	700	1.6	23.0	Y = 0.984x + 0.113	0.998

Point-of-care studies

Point-of-care studies have been performed at three individual district health centers in primary care with both nursing and physician services. The studies were performed on capillary samples in duplicates tested on HemoCue Hb 301 vs HemoCue Hb 201⁺ (1). The Hb 201⁺ was verified against the International Council for Standardization in Haematology method (ICSH)³ using venous blood (2).

Study	N	Min g/dL	Max g/dL	Regression line	Correlation coefficient (r)
1	186	8.3	18.0	$Y = 0.897x + 1.357$	0.930
2	700	1.6	23.0	$Y = 1.005x - 0.041$	0.999

Technical Specifications

Dimensions: 160 x 140 x 70 mm (6.29 x 5.51 x 2.76 inches)

Weight: 500 g (1.10 pounds) (with batteries installed)

4 batteries type AA

Power adapter: CE marked

Only use adapters, as listed under "Adapters".

Pollution degree: 2

Overvoltage category: II

Atmospheric pressure: 700 hPa to 1060 hPa.

Equipment not suitable for use in the presence of flammable mixtures.

The instrument is tested according to

- IEC/EN 61010-1, 2001
- IEC/EN 60601-1 2:nd edition clause 19, 20 and 57.10 (Leakage current, dielectric strength and insulation distance)
- IEC/EN 60601-1-2
- IEC/EN 61326

and complies with the IVD Medical Device Directive 98/79/EC.

The instrument is made for continuous mode.

Adapters

Country: EU/US/GB

Type: FW7333SM/12

Input: 100 V~-240 V~/50-60 Hz/<200 mA

Country: EU/US/GB

Type: HCA01

Input: 100 V~-240 V~/50-60 Hz/<500 mA

Warnings

The device is tested according to standard IEC/EN 60601-1-2 and is found to comply with the standard.

Despite this compliance it is impossible to predict any possible effects by other nearby standing instruments, (stationary, portable or mobile units) or the possible impact of electromagnetic radiance. This is the reason we need to inform users of this analyzer that disturbance from other equipment might affect the performance of the analyzer. Should you note that this is the fact, please contact your local HemoCue distributor.

The HemoCue Hb 301 is intended for use in the electromagnetic environment specified in Technical specifications. The customer or user of the HemoCue Hb 301 system should assure that it is used in such an environment. The HemoCue Hb 301 Analyzer uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.

The HemoCue Hb 301 system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (i.e. IEC 60950-1 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all configurations must comply with the system standard, IEC 60601-1-1. Connections of additional equipment to the signal input or signal output connections is deemed to be "configuring a medical system", and therefore assumes responsibility

for the systems compliance within the requirements of IEC 60601-1-1. For further information, please see the attachment Technical Specifications (EMC-RF) or contact your local HemoCue distributor.

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue® Hb 301 Analyzer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds For explanation of U see NOTE 1	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> <p>See NOTE 2 and NOTE 3</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d=1.2\sqrt{P}$</p> <p>$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 KHz to 80 Mhz, field strength should be less than 3 V/m.

Technical specifications (EMC-RF)**Use only cables with the following specification:**

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

Guidance and manufacturer's declaration – electromagnetic emissions		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Warranty

The analyzer carries a 24-month warranty from the day of receipt. After the expiration date of the warranty, maintenance and repairs are offered at a fixed price.

Service and disposal

The analyzer should be cleaned as recommended under "Maintenance" prior to service or disposal. Consult local environmental authorities for proper disposal.

Spare parts and accessories:

The following accessories and spare parts are available:

- Power adapter
- Battery lid
- Cuvette holder
- HemoCue Cleaners
- HemoCue Lancets

Patent

This product is covered by US patent no. 6 831 733.

Symbols used

Attention, see instructions for use

CE mark

Class II equipment

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

Serial port

Relative humidity, non-condensing

References

1. HemoCue Hb 301 Microcuvettes package insert¹.
2. Dacie and Lewis, Practical Haematology Ninth edition.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.
5. NCCLS Document EP6-A Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.

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Humedad relativa, no condensada

Referencias

1. respecto del producto HemoCue Hb 301 Microcuvettes¹.
2. Dacie and Lewis, Practical Haematology Ninth edition.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.
5. NCCLS Document EP6-A Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.

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Porta serial

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Referências

1. Bula das HemoCue Hb 301 Microcuvettes¹.
2. Dacie and Lewis, Practical Haematology Ninth edition.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.
5. NCCLS Document EP6-A Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.

Fabricante

HemoCue AB

Box 1204

SE-262 23 Ängelholm

Suécia

Telephone: +46 77 570 0210

Fax: +46 77 570 0212

E-mail: info@hemocue.se

www.hemocue.com

HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 77 570 02 10. Fax +46 77 570 02 12.
info@hemocue.se www.hemocue.com